

Mouse Bcl10/B-cell lymphoma/leukemia 10 ELISA Kit

Catalog No: #EK20079

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Bcl10/B-cell lymphoma/leukemia 10 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse B-cell lymphoma/leukemia 10
Target Name	B-cell lymphoma/leukemia 10
Other Names	B-cell CLL/lymphoma 10, Bcl-10, CARD-containing molecule enhancing NF-kappa-B, CARD-like apoptotic protein, mCLAP, CED-3/ICH-1 prodomain homologous E10-like regulator, mCIPER, Cellular homolog of vCARMEN, cCARMEN, Cellular-E10, c-E10, Mammalian CARD-containing adapter molecule E10, mE10, Ciper, Clap
Accession No.	Q9Z0H7
Uniprot	Q9Z0H7
GeneID	12042;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.174ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bcl10

Note: This product is for in vitro research use only