

Rat Bmf/Bcl-2-modifying factor ELISA Kit

Catalog No: #EK20081

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Bmf/Bcl-2-modifying factor ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Bcl-2-modifying factor
Target Name	Bcl-2-modifying factor
Accession No.	Q8K589
Uniprot	Q8K589
GeneID	246142;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.175ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Bmf

Note: This product is for in vitro research use only