

Mouse Brms1/Breast cancer metastasis-suppressor 1 homolog ELISA Kit



Catalog No: #EK20083

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Mouse Brms1/Breast cancer metastasis-suppressor 1 homolog ELISA Kit |
| Brief Description | ELISA Kits |
| Applications | ELISA |
| Species Reactivity | Ms |
| Specificity | Natural and recombinant Mouse Breast cancer metastasis-suppressor 1 homolog |
| Target Name | Breast cancer metastasis-suppressor 1 homolog |
| Accession No. | Q99N20 |
| Uniprot | Q99N20 |
| GeneID | 107392; |

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Brms1

Note: This product is for in vitro research use only