

Rat Brms1/Breast cancer metastasis-suppressor 1 homolog ELISA Kit

Catalog No: #EK20084

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Brms1/Breast cancer metastasis-suppressor 1 homolog ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Breast cancer metastasis-suppressor 1 homolog
Target Name	Breast cancer metastasis-suppressor 1 homolog
Accession No.	Q5M7T3
Uniprot	Q5M7T3
GeneID	108348098;293668;

Application Details

Assay Type: SandwichDetect Range:31.2-2000pg/mL Sensitivity:15.69pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nmgene_full_name:Brms1

Note: This product is for in vitro research use only