

Human CASP6/Caspase-6 ELISA Kit

Catalog No: #EK20085

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CASP6/Caspase-6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Caspase-6
Target Name	Caspase-6
Other Names	CASP-6, Apoptotic protease Mch-2, MCH2
Accession No.	P55212
Uniprot	P55212
GeneID	839;

Application Details

Assay Type: Sandwich Detect Range: 31.2-2000pg/mL Sensitivity: 15.71pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: CASP6

Note: This product is for in vitro research use only