

Mouse Casp6/Caspase-6 ELISA Kit

Catalog No: #EK20086

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Casp6/Caspase-6 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Caspase-6
Target Name	Caspase-6
Other Names	CASP-6, Apoptotic protease Mch-2
Accession No.	O08738
Uniprot	O08738
GeneID	12368;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.156ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Casp6

Note: This product is for in vitro research use only