

Mouse Cda/Cytidine deaminase ELISA Kit

Catalog No: #EK20089

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cda/Cytidine deaminase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cytidine deaminase
Target Name	Cytidine deaminase
Other Names	Cytidine aminohydrolase,Cdd
Accession No.	P56389
Uniprot	P56389
GeneID	72269;

Application Details

Assay Type: Sandwich Detect Range:78-5000pg/ml Sensitivity:39.8pg/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Cda

Note: This product is for in vitro research use only