

Mouse Cldn1/Claudin-1 ELISA Kit

Catalog No: #EK20090

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cldn1/Claudin-1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Claudin-1
Target Name	Claudin-1
Accession No.	O88551
Uniprot	O88551
GeneID	12737;

Application Details

Assay Type: Sandwich Detect Range: 0.156-10ng/ml Sensitivity: 0.078ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cldn1

Note: This product is for in vitro research use only