

Human CRK/Adapter molecule crk ELISA Kit

Catalog No: #EK20094

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CRK/Adapter molecule crk ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Adapter molecule crk
Target Name	Adapter molecule crk
Other Names	Proto-oncogene c-Crk,p38
Accession No.	P46108
Uniprot	P46108
GeneID	1398;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.078ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CRK

Note: This product is for in vitro research use only