

Rat Crk/Adapter molecule crk ELISA Kit

Catalog No: #EK20095

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Rat Crk/Adapter molecule crk ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Rt
Specificity	Natural and recombinant Rat Adapter molecule crk
Target Name	Adapter molecule crk
Other Names	Proto-oncogene c-Crk,p38,Crko
Accession No.	Q63768
Uniprot	Q63768
GeneID	54245;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.097 ng/ml Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:Crk

Note: This product is for in vitro research use only