

Bovine CRIP1/Cysteine-rich protein 1 ELISA Kit

Catalog No: #EK20096

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bovine CRIP1/Cysteine-rich protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Bovine
Specificity	Natural and recombinant Bovine Cysteine-rich protein 1
Target Name	Cysteine-rich protein 1
Other Names	CRP-1
Accession No.	Q56K04
Uniprot	Q56K04
GeneID	574093;

Application Details

Assay Type: Sandwich Detect Range:0.625-40ng/mL Sensitivity:0.308ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CRIP1

Note: This product is for in vitro research use only