

Human CRIP1/Cysteine-rich protein 1 ELISA Kit

Catalog No: #EK20097

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CRIP1/Cysteine-rich protein 1 ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cysteine-rich protein 1
Target Name	Cysteine-rich protein 1
Other Names	CRP-1,Cysteine-rich heart protein,CRHP,Cysteine-rich intestinal protein,CRIP,CRIP,CRP1
Accession No.	P50238
Uniprot	P50238
GeneID	1396;

Application Details

Assay Type: Sandwich Detect Range:0.156-10ng/mL Sensitivity:0.082ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CRIP1

Note: This product is for in vitro research use only