

## Human CTSE/Cathepsin E ELISA Kit

Catalog No: #EK20099

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CTSE/Cathepsin E ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Cathepsin E
Target Name	Cathepsin E
Accession No.	P14091
Uniprot	P14091
GeneID	1510;

## Application Details

Assay Type: Sandwich Detect Range: 78-5000pg/ml Sensitivity: 40.6pg/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.  
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene\_full\_name: CTSE

Note: This product is for in vitro research use only