

Mouse Cubn/Cubilin ELISA Kit

Catalog No: #EK20100

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mouse Cubn/Cubilin ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Ms
Specificity	Natural and recombinant Mouse Cubilin
Target Name	Cubilin
Other Names	Intrinsic factor-cobalamin receptor, Ifcr
Accession No.	Q9JLB4
Uniprot	Q9JLB4
GeneID	65969;

Application Details

Assay Type: Sandwich Detect Range: 0.312-20ng/mL Sensitivity: 0.11ng/mL Sample Type: Serum, plasma, tissue homogenates and other biological fluids.
Sample Volume: 100ul Assay Time: 3-5h Detection wavelength: 450 nm gene_full_name: Cubn

Note: This product is for in vitro research use only