

Human CRAT/Carnitine O-acetyltransferase ELISA Kit

Catalog No: #EK20101

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human CRAT/Carnitine O-acetyltransferase ELISA Kit
Brief Description	ELISA Kits
Applications	ELISA
Species Reactivity	Hu
Specificity	Natural and recombinant Human Carnitine O-acetyltransferase
Target Name	Carnitine O-acetyltransferase
Other Names	Carnitine acetylase,Carnitine acetyltransferase,CAT,CAT1
Accession No.	P43155
Uniprot	P43155
GeneID	1384;

Application Details

Assay Type: Sandwich Detect Range:0.312-20ng/mL Sensitivity:0.156ng/mL Sample Type:Serum, plasma, tissue homogenates and other biological fluids.

Sample Volume:100ul Assay Time:3-5h Detection wavelength:450 nm gene_full_name:CRAT

Note: This product is for in vitro research use only