

LPPL Antibody

Catalog No: #44728

Package Size: #44728-1 50ul #44728-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	LPPL Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	LPPL antibody detects endogenous levels of total LPPL
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human LPPL
Target Name	LPPL
Other Names	CLC
Accession No.	Swiss-Prot#: Q05315 NCBI Gene ID:1178
Uniprot	Q05315
GeneID	1178;
Calculated MW	16 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only