

ISK7 Antibody

Catalog No: #44731

Package Size: #44731-1 50ul #44731-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ISK7 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	ISK7antibody detects endogenous levels of total?ISK7
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ISK7
Target Name	ISK7
Other Names	Serine protease inhibitor Kazal-type 7,ECRG-2,ECG2,SPINK7
Accession No.	Swiss-Prot#: P58062NCBI Gene ID:84651
Uniprot	P58062
GeneID	84651;
Calculated MW	9 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Images



Western Blot analysis of lysates of HeLa cell and 293T cell, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only