

MRCKG Antibody

Catalog No: #44745

Package Size: #44745-1 50ul #44745-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	MRCKG Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	MRCKG antibody detects endogenous levels of total MRCKG
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MRCKG
Target Name	MRCKG
Other Names	CDC42BPG
Accession No.	Swiss-Prot#: Q6DT37 NCBI Gene ID:55561
Uniprot	Q6DT37
GeneID	55561;
Calculated MW	20 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only