

## ZNF497 Antibody

Catalog No: #44746

Package Size: #44746-1 50ul #44746-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ZNF497 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	ZNF497 antibody detects endogenous levels of total ZNF497
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human ZNF497
Target Name	ZNF497
Other Names	Zinc finger protein 497
Accession No.	Swiss-Prot#: Q6ZNH5 NCBI Gene ID:162968
Uniprot	Q6ZNH5
GeneID	162968;
Calculated MW	55 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

## Images



Western Blot analysis of HeLa lysates, using primary antibody at 1:1000 dilution.

---

Note: This product is for in vitro research use only