

## NUFIP2 Antibody

Catalog No: #44748

Package Size: #44748-1 50ul #44748-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	NUFIP2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	NUFIP2 antibody detects endogenous levels of total NUFIP2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human NUFIP2
Target Name	NUFIP2
Other Names	82 kDa FMRP-interacting protein antibody 82-FIP antibody Cell proliferation-inducing gene 1 protein antibody FMRP-interacting protein 2 antibody Nuclear fragile X mental retardation-interacting protein 2 antibody NUFIP2 antibody NUFP2_HUMAN antibody
Accession No.	Swiss-Prot#: Q7Z417 NCBI Gene ID:57532
Uniprot	Q7Z417
GeneID	57532;
Calculated MW	76 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

## Images



Western Blot analysis of HeLa lysates, using primary antibody at 1:1000 dilution.

---

Note: This product is for in vitro research use only