

HP1 α (Phospho-Ser92) Conjugated Antibody

Catalog No: #C11777



Package Size: #C11777-AF350 100ul #C11777-AF405 100ul #C11777-AF488 100ul
 #C11777-AF555 100ul #C11777-AF594 100ul #C11777-AF647 100ul
 #C11777-AF680 100ul #C11777-AF750 100ul #C11777-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	HP1 α (Phospho-Ser92) Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of HP1 α only when phosphorylated at serine 92.
Immunogen Description	Peptide sequence around phosphorylation site of Serine 92(R-S-S(p)-N-F) derived from Human HP1 α .
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	HP1A; HP1 alpha;antigen p25
Accession No.	Swiss-Prot#:P45973NCBI Gene ID:23468NCBI mRNA#:NM_001127321.1. NCBI Protein#:NP_001120793.1.
Uniprot	P45973
GeneID	23468;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	22
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.

Background

Heterochromatin protein-1 (HP1) is a methyl-lysine binding protein localized at heterochromatin sites, where it mediates gene silencing.

Note: This product is for in vitro research use only