

DIRC1 Antibody

Catalog No: #44762

Package Size: #44762-1 50ul #44762-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | DIRC1 Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Immunogen affinity purified |
| Applications | WB |
| Species Reactivity | Hu |
| Specificity | DIRC1 antibody detects endogenous levels of total DIRC1 |
| Immunogen Type | Peptide |
| Immunogen Description | A synthesized peptide derived from human DIRC1 |
| Target Name | DIRC1 |
| Other Names | Disrupted in renal carcinoma protein 1, Disrupted in renal cancer protein |
| Accession No. | Swiss-Prot#: Q969H9 NCBI Gene ID: 116093 |
| Uniprot | Q969H9 |
| GeneID | 116093; |
| Calculated MW | 11 kDa |
| Concentration | 1.0mg/mL |
| Formulation | Liquid in PBS containing 50% glycerol, and 0.02% sodium azide. |
| Storage | Store at -20°C |

Application Details

Western blotting: □ 1:500-1:2000

Images



Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only