

DMRTB1 Antibody

Catalog No: #44764

Package Size: #44764-1 50ul #44764-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DMRTB1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	DMRTB1 antibody detects endogenous levels of total DMRTB1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human DMRTB1
Target Name	DMRTB1
Other Names	DMRT B 1 antibody DMRT like family B with proline rich C terminal 1 antibody DMRTB_HUMAN antibody Dmrtb1 antibody Doublesex- and mab-3-related transcription factor B1 antibody FLJ32711 antibody FLJ40434 antibody HGNC:13913 antibody
Accession No.	Swiss-Prot#: Q96MA1 NCBI Gene ID:63948
Uniprot	Q96MA1
GeneID	63948;
Calculated MW	44 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only