

AP-1/Jun D (Phospho-Ser73/100) Polyclonal Conjugated Antibody

Catalog No: #C12181

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C12181-AF350 100ul #C12181-AF405 100ul #C12181-AF488 100ul

#C12181-AF555 100ul #C12181-AF594 100ul #C12181-AF647 100ul

#C12181-AF680 100ul #C12181-AF750 100ul #C12181-Biotin 100ul

Description

Product Name	AP-1/Jun D (Phospho-Ser73/100) Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	Phospho-AP-1 (S73) Polyclonal Antibody detects endogenous levels of AP-1 protein only when phosphorylated at S73.
Immunogen Description	Synthesized peptide derived from human AP-1 around the phosphorylation site of S73.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	JUN; Transcription factor AP-1; Activator protein 1; AP1; Proto-oncogene c-Jun; V-jun avian sarcoma virus 17 oncogene homolog; p39; JUND; Transcription factor jun-D
Accession No.	Swiss-Prot#:P05412/P17535NCBI Gene ID:3725/3728
Uniprot	P05412
GeneID	3725;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	36
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only