

CREB-2 (Phospho-Ser219) Polyclonal Conjugated Antibody



Catalog No: #C12265

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C12265-AF350 100ul #C12265-AF405 100ul #C12265-AF488 100ul

#C12265-AF555 100ul #C12265-AF594 100ul #C12265-AF647 100ul

#C12265-AF680 100ul #C12265-AF750 100ul #C12265-Biotin 100ul

Description

Product Name	CREB-2 (Phospho-Ser219) Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	Phospho-CREB-2 (S219) Polyclonal Antibody detects endogenous levels of CREB-2 protein only when phosphorylated at S219.
Immunogen Description	Synthesized peptide derived from human CREB-2 around the phosphorylation site of S219.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ATF4; CREB2; TXREB; Cyclic AMP-dependent transcription factor ATF-4; cAMP-dependent transcription factor ATF-4; Activating transcription factor 4; Cyclic AMP-responsive element-binding protein 2; CREB-2; cAMP-responsive element-binding prot
Accession No.	Swiss-Prot#:P18848NCBI Gene ID:468
Uniprot	P18848
GeneID	468;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	38
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only