

# HSP90 $\beta$ (Phospho-Ser226) Polyclonal Conjugated Antibody

Catalog No: #C12277

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #C12277-AF350 100ul #C12277-AF405 100ul #C12277-AF488 100ul

#C12277-AF555 100ul #C12277-AF594 100ul #C12277-AF647 100ul

#C12277-AF680 100ul #C12277-AF750 100ul #C12277-Biotin 100ul

## Description

Product Name	HSP90 $\beta$ (Phospho-Ser226) Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	Phospho-HSP90 $\beta$ (S226) Polyclonal Antibody detects endogenous levels of HSP90 $\beta$ protein only when phosphorylated at S226.
Immunogen Description	Synthesized peptide derived from human HSP90 $\beta$ around the phosphorylation site of S226.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	HSP90AB1; HSP90B; HSPC2; HSPCB; Heat shock protein HSP 90-beta; HSP 90; Heat shock 84 kDa; HSP 84; HSP84
Accession No.	Swiss-Prot#:P08238NCBI Gene ID:3326
Uniprot	P08238
GeneID	3326;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	83
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

---

Note: This product is for in vitro research use only