

MEK Kinase-1 (Phospho-Thr1402) Polyclonal Conjugated Antibody

Catalog No: #C12336

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C12336-AF350 100ul #C12336-AF405 100ul #C12336-AF488 100ul

#C12336-AF555 100ul #C12336-AF594 100ul #C12336-AF647 100ul

#C12336-AF680 100ul #C12336-AF750 100ul #C12336-Biotin 100ul

Description

Product Name	MEK Kinase-1 (Phospho-Thr1402) Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	Phospho-MEK Kinase-1 (T1402) Polyclonal Antibody detects endogenous levels of MEK Kinase-1 protein only when phosphorylated at T1402.
Immunogen Description	Synthesized peptide derived from human MEK Kinase-1 around the phosphorylation site of T1402.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	MAP3K1; MAPKKK1; MEKK; MEKK1; Mitogen-activated protein kinase kinase kinase 1; MAPK/ERK kinase kinase 1; MEK kinase 1; MEKK 1
Accession No.	Swiss-Prot#:Q13233NCBI Gene ID:4214
Uniprot	Q13233
GeneID	4214;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	130
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only