

Gab 1 (Phospho-Tyr659) Polyclonal Conjugated Antibody

Catalog No: #C12359

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C12359-AF350 100ul #C12359-AF405 100ul #C12359-AF488 100ul

#C12359-AF555 100ul #C12359-AF594 100ul #C12359-AF647 100ul

#C12359-AF680 100ul #C12359-AF750 100ul #C12359-Biotin 100ul

Description

Product Name	Gab 1 (Phospho-Tyr659) Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	Phospho-Gab 1 (Y659) Polyclonal Antibody detects endogenous levels of Gab 1 protein only when phosphorylated at Y659.
Immunogen Description	Synthesized peptide derived from human Gab 1 around the phosphorylation site of Y659.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	GAB1; GRB2-associated-binding protein 1; GRB2-associated binder 1; Growth factor receptor bound protein 2-associated protein 1
Accession No.	Swiss-Prot#:Q13480NCBI Gene ID:2549
Uniprot	Q13480
GeneID	2549;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	76
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only