

PDGFR- β (Phospho-Tyr1021) Polyclonal Conjugated Antibody



Catalog No: #C12374

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C12374-AF350 100ul #C12374-AF405 100ul #C12374-AF488 100ul

#C12374-AF555 100ul #C12374-AF594 100ul #C12374-AF647 100ul

#C12374-AF680 100ul #C12374-AF750 100ul #C12374-Biotin 100ul

Description

Product Name	PDGFR- β (Phospho-Tyr1021) Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	Phospho-PDGFR- β (Y1021) Polyclonal Antibody detects endogenous levels of PDGFR- β protein only when phosphorylated at Y1021.
Immunogen Description	Synthesized peptide derived from human PDGFR- β around the phosphorylation site of Y1021.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	PDGFRB; PDGFR; PDGFR1; Platelet-derived growth factor receptor beta; PDGF-R-beta; PDGFR-beta; Beta platelet-derived growth factor receptor; Beta-type platelet-derived growth factor receptor; CD140 antigen-like family member B; Platelet-deri
Accession No.	Swiss-Prot#:P09619NCBI Gene ID:5159
Uniprot	P09619
GeneID	5159;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only