

## TOPRS Antibody

Catalog No: #44788

Package Size: #44788-1 50ul #44788-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	TOPRS Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	TOPRSantibody detects endogenous levels of total?TOPRS
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TOPRS
Target Name	TOPRS
Other Names	LUN, TP53BPL, TOPORS
Accession No.	Swiss-Prot#: Q9NS56NCBI Gene ID:10210
Uniprot	Q9NS56
GeneID	10210;
Calculated MW	119 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

## Images



Western Blot analysis of HeLa lysates, using primary antibody at 1:1000 dilution.

---

Note: This product is for in vitro research use only