

Sodium Potassium ATPase alpha-1 (Phospho-Tyr260) Conjugated Antibody

Catalog No: #C12398

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C12398-AF350 100ul #C12398-AF405 100ul #C12398-AF488 100ul

#C12398-AF555 100ul #C12398-AF594 100ul #C12398-AF647 100ul

#C12398-AF680 100ul #C12398-AF750 100ul #C12398-Biotin 100ul

Description

Product Name	Sodium Potassium ATPase alpha-1 (Phospho-Tyr260) Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	The antibody detects endogenous Sodium Potassium ATPase alpha-1 when Phospho occurs at Tyr260
Immunogen Description	Synthetic peptide from human protein at AA range: 230-290
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Sodium/potassium-transporting ATPase subunit alpha-1 (Na(+)/K(+) ATPase alpha-1 subunit) (EC 3.6.3.9) (Sodium pump subunit alpha-1)
Accession No.	Swiss-Prot#:P05023NCBI Gene ID:476
Uniprot	P05023
GeneID	476;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	115
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only