

## TSSC4 Antibody

Catalog No: #44796

Package Size: #44796-1 50ul #44796-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TSSC4 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	TSSC4antibody detects endogenous levels of total TSSC4
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TSSC4
Target Name	TSSC4
Other Names	TSSC4
Accession No.	Swiss-Prot#: Q9Y5U2NCBI Gene ID:10078
Uniprot	Q9Y5U2
GeneID	10078;
Calculated MW	34 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

## Images



---

Note: This product is for in vitro research use only