

c-Jun(Phospho-T91) Conjugated Antibody

Catalog No: #C13420



Package Size: #C13420-AF350 100ul #C13420-AF405 100ul #C13420-AF488 100ul
 #C13420-AF555 100ul #C13420-AF594 100ul #C13420-AF647 100ul
 #C13420-AF680 100ul #C13420-AF750 100ul #C13420-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	c-Jun(Phospho-T91) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Species Reactivity	Hu
Immunogen Description	recombinant protein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Ah receptor antibody AhR antibody AHR_HUMAN antibody Aromatic hydrocarbon receptor antibody Aryl hydrocarbon receptor antibody Aryl hydrocarbon receptor precursor antibody bHLHe76 antibody Class E basic helix loop helix protein 76 antibody Class E basic helix-loop-helix protein 76 antibody HGNC:348 antibody
Accession No.	Swiss-Prot#:P35869
Uniprot	P35869
GeneID	196;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	96
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250
 AF405 conjugated: most applications: 1: 50 - 1: 250
 AF488 conjugated: most applications: 1: 50 - 1: 250
 AF555 conjugated: most applications: 1: 50 - 1: 250
 AF594 conjugated: most applications: 1: 50 - 1: 250
 AF647 conjugated: most applications: 1: 50 - 1: 250
 AF680 conjugated: most applications: 1: 50 - 1: 250
 AF750 conjugated: most applications: 1: 50 - 1: 250

Background

2,3,7,8-tetrachlorodibenzo-p-dioxin (TCDD) is the prototype for a family of toxic halogenated aromatic compounds that are thought to cause adverse reproductive, immunologic and metabolic effects. Many biological responses to TCDD are mediated through ligand binding to the aromatic hydrocarbon (Ah) receptor, also known as AhR. Ah Receptor is a ligand dependent transcription factor that interacts with specific DNA sequences, termed xenobiotic responsive elements (XREs), and that lies upstream of TCDD-inducible genes. Upon binding to the ligand, Ah Receptor binds to the Ah Receptor nuclear translocator (Arnt), and the complex is translocated from the cytoplasm to the nucleus. Arnt is required for Ah Receptor to bind to XRE. Ah Receptor and Arnt are members of a family of transcription factors that contain a basic helix-loop-helix motif and a common "PAS" motif.

Note: This product is for in vitro research use only