

PLK1(Phospho-T210) Conjugated Antibody

Catalog No: #C13421



Package Size: #C13421-AF350 100ul #C13421-AF405 100ul #C13421-AF488 100ul
 #C13421-AF555 100ul #C13421-AF594 100ul #C13421-AF647 100ul
 #C13421-AF680 100ul #C13421-AF750 100ul #C13421-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	PLK1(Phospho-T210) Conjugated Antibody
Host Species	Rabbit
Clonality	Monoclonal
Immunogen Description	recombinant protein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	HA epitope tag antibody HA1 antibody HA2 antibody hemagglutinin antibody Hemagglutinin HA1 chain antibody Hemagglutinin HA2 chain antibody
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

Human influenza hemagglutinin (HA) is a surface glycoprotein required for the infectivity of the human virus. The HA tag is derived from the HA molecule corresponding to amino acids 98-106. This antibody is used to detect proteins that are tagged with HA epitope (YPYDVPDYA) and are expressed in prokaryotic and eukaryotic cells.

Note: This product is for in vitro research use only