

## HSPB3 Antibody

Catalog No: #44808

Package Size: #44808-1 50ul #44808-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	HSPB3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	HSPB3antibody detects endogenous levels of total?HSPB3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human HSPB3
Target Name	HSPB3
Other Names	2310035K17Rik antibody A1844863 antibody DHMN2C antibody Heat shock 17 kDa protein antibody Heat shock 27kD protein family, member 3; antibody Heat shock 27kDa protein 3 antibody Heat shock protein 3 antibody Heat shock protein beta 3 antibody HMN2C antibody HSP 17 antibody HSPL27 antibody Protein 3
Accession No.	Swiss-Prot#: Q12988NCBI Gene ID:8988
Uniprot	Q12988
GeneID	8988;
Calculated MW	17 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

## Application Details

Western blotting: □ 1:500-1:2000

## Images



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only