

Akt2 (Ab-474) Conjugated Antibody

Catalog No: #C21155



Package Size: #C21155-AF350 100ul #C21155-AF405 100ul #C21155-AF488 100ul
 #C21155-AF555 100ul #C21155-AF594 100ul #C21155-AF647 100ul
 #C21155-AF680 100ul #C21155-AF750 100ul #C21155-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Akt2 (Ab-474) Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu Ms Rt
Specificity	The antibody detects endogenous level of total Akt2 protein.
Immunogen Description	Peptide sequence around aa. 472~476 (Q-F-S-Y-S) derived from Human Akt2.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	PKB beta;Akt-2; Protein kinase B;RAC-PK-beta
Accession No.	Swiss-Prot#:P31751NCBI Gene ID:208NCBI mRNA#:NM_001626.3 NCBI Protein#: NP_001617.1
Uniprot	P31751
GeneID	208;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	60
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Product Description

Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.

Background

General protein kinase capable of phosphorylating several known proteins.

Note: This product is for in vitro research use only