

IGFL1 Antibody

Catalog No: #44813

Package Size: #44813-1 50ul #44813-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	IGFL1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	IGFL1 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human IGFL1
Target Name	IGFL1
Other Names	APRG644 antibody IGFL1 antibody IGFL1_HUMAN antibody Insulin growth factor-like family member 1 antibody UNQ644 antibody
Accession No.	Swiss-Prot#: Q6UW32NCBI Gene ID:374918
Uniprot	Q6UW32
GeneID	374918;
Calculated MW	12 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only