

## RE57611p Conjugated Antibody

Catalog No: #C21416



Package Size: #C21416-AF350 100ul #C21416-AF405 100ul #C21416-AF488 100ul  
 #C21416-AF555 100ul #C21416-AF594 100ul #C21416-AF647 100ul  
 #C21416-AF680 100ul #C21416-AF750 100ul #C21416-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	RE57611p Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Fruit fly
Specificity	The antibody detects endogenous level of total RE57611 protein.
Immunogen Description	Peptide sequence around aa.178~182 (K-R-Q-T-L) derived from Fruit fly RE57611p.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Accession No.	Swiss-Prot#:D3DMZ6 NCBI Protein#:ADC54212
Uniprot	D3DMZ6
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	23
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

## Product Description

---

Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.

## Background

---

CG8889 gene product from transcript CG8889-RA

---

Note: This product is for in vitro research use only