

# Protein Phosphatase 4C Conjugated Monoclonal Antibody

Catalog No: #C27020

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #C27020-AF350 100ul #C27020-AF405 100ul #C27020-AF488 100ul

#C27020-AF555 100ul #C27020-AF594 100ul #C27020-AF647 100ul

#C27020-AF680 100ul #C27020-AF750 100ul #C27020-Biotin 100ul

## Description

Product Name	Protein Phosphatase 4C Conjugated Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Immunogen Description	Recombinant full length protein (307aa) of human PPP4C.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	serine/threonine-protein phosphatase 4 catalytic subunit PP-X protein phosphatase X, catalytic subunit protein phosphatase 4 (formerly X), catalytic subunit
Accession No.	Swiss-Prot#:P60510 NCBI Gene ID:5531
Uniprot	P60510
GeneID	5531;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	34
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

---

Note: This product is for in vitro research use only