

## Keratin 7(C-term) Conjugated Monoclonal Antibody

Catalog No: #C27028



Package Size: #C27028-AF350 100ul #C27028-AF405 100ul #C27028-AF488 100ul  
 #C27028-AF555 100ul #C27028-AF594 100ul #C27028-AF647 100ul  
 #C27028-AF680 100ul #C27028-AF750 100ul #C27028-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Keratin 7(C-term) Conjugated Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Immunogen Description	Recombinant fragment, , corresponding to amino acids 256aa-415aa of human KRT7 (C-terminus).
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	keratin, type II cytoskeletal 7 CK-7 keratin-7 sarcolectin cytokeratin 7 cytokeratin-7 type-II keratin Kb7 type II mesothelial keratin K7 keratin, 55K type II cytoskeletal keratin, simple epithelial type I, K7
Accession No.	Swiss-Prot#:P08729 NCBI Gene ID:3855
Uniprot	P08729
GeneID	3855;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	55
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

Note: This product is for in vitro research use only