

TRIM33 Conjugated Monoclonal Antibody

Catalog No: #C27032



Package Size: #C27032-AF350 100ul #C27032-AF405 100ul #C27032-AF488 100ul
 #C27032-AF555 100ul #C27032-AF594 100ul #C27032-AF647 100ul
 #C27032-AF680 100ul #C27032-AF750 100ul #C27032-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	TRIM33 Conjugated Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Immunogen Description	Recombinant fragment, corresponding to amino acids 712aa-853aa of human TRIM33 (C-terminus).
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	E3 ubiquitin-protein ligase TRIM33 TIF1-gamma protein Rfg7 ectodermin homolog RET-fused gene 7 protein tripartite motif-containing 33 transcriptional intermediary factor 1 gamma
Accession No.	Swiss-Prot#:Q9UPN9 NCBI Gene ID:51592
Uniprot	Q9UPN9
GeneID	51592;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	140
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only