

TBLR1 Conjugated Monoclonal Antibody

Catalog No: #C27036



Package Size: #C27036-AF350 100ul #C27036-AF405 100ul #C27036-AF488 100ul
 #C27036-AF555 100ul #C27036-AF594 100ul #C27036-AF647 100ul
 #C27036-AF680 100ul #C27036-AF750 100ul #C27036-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	TBLR1 Conjugated Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Immunogen Description	Recombinant full length protein (514aa) of human TBLR1.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	F-box-like/WD repeat-containing protein TBL1XR1 TBL1-related protein 1 transducin beta-like 1X-related protein 1 nuclear receptor co-repressor/HDAC3 complex subunit nuclear receptor co-repressor/HDAC3 complex subunit TBLR1
Accession No.	Swiss-Prot#:Q9BZK7 NCBI Gene ID:79718
Uniprot	Q9BZK7
GeneID	79718;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	60
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250
 AF405 conjugated: most applications: 1: 50 - 1: 250
 AF488 conjugated: most applications: 1: 50 - 1: 250
 AF555 conjugated: most applications: 1: 50 - 1: 250
 AF594 conjugated: most applications: 1: 50 - 1: 250
 AF647 conjugated: most applications: 1: 50 - 1: 250
 AF680 conjugated: most applications: 1: 50 - 1: 250
 AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only