

Ikaros(C-term) Conjugated Monoclonal Antibody

Catalog No: #C27062



Package Size: #C27062-AF350 100ul #C27062-AF405 100ul #C27062-AF488 100ul
 #C27062-AF555 100ul #C27062-AF594 100ul #C27062-AF647 100ul
 #C27062-AF680 100ul #C27062-AF750 100ul #C27062-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ikaros(C-term) Conjugated Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Immunogen Description	Recombinant fragment, corresponding to amino acids 249-415aa of human IKAROS(C-terminus).
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	DNA-binding protein Ikaros CLL-associated antigen KW-6 Ikaros (zinc finger protein) ikaros family zinc finger protein 1 lymphoid transcription factor LyF-1 zinc finger protein, subfamily 1A, 1 (Ikaros
Accession No.	Swiss-Prot#:Q13422 NCBI Gene ID:10320
Uniprot	Q13422
GeneID	10320;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	58
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only