

Dynactin 1(N-terminus) Conjugated Monoclonal Antibody

Catalog No: #C27082

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C27082-AF350 100ul #C27082-AF405 100ul #C27082-AF488 100ul

#C27082-AF555 100ul #C27082-AF594 100ul #C27082-AF647 100ul

#C27082-AF680 100ul #C27082-AF750 100ul #C27082-Biotin 100ul

Description

Product Name	Dynactin 1(N-terminus) Conjugated Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Immunogen Description	Purified recombinant human DCTN1 protein fragments expressed in E.coli.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	P135, DP-150, DAP-150
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	150
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

Required for the cytoplasmic dynein-driven retrograde movement of vesicles and organelles along microtubules. Dynein-dynactin interaction is a key

component of the mechanism of axonal transport of vesicles and organelles.

Note: This product is for in vitro research use only