

CNNM1 Antibody

Catalog No: #44845

Package Size: #44845-1 50ul #44845-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CNNM1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	CNNM1 antibody detects endogenous levels of total CNNM1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CNNM1
Target Name	CNNM1
Other Names	ACDP1 antibody Ancient conserved domain protein 1 antibody Ancient conserved domain-containing protein 1 antibody CNNM1 antibody CNNM1_HUMAN antibody Cyclin-M1 antibody FLJ31632 antibody Metal transporter CNNM1 antibody OTTHUMP00000020253 antibody
Accession No.	Swiss-Prot#: Q9NRU3 NCBI Gene ID:26507
Uniprot	Q9NRU3
GeneID	26507;
Calculated MW	72kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Note: This product is for in vitro research use only