

## UCHL1/PGP9.5 Conjugated Monoclonal Antibody

Catalog No: #C27188



Package Size: #C27188-AF350 100ul #C27188-AF405 100ul #C27188-AF488 100ul  
 #C27188-AF555 100ul #C27188-AF594 100ul #C27188-AF647 100ul  
 #C27188-AF680 100ul #C27188-AF750 100ul #C27188-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	UCHL1/PGP9.5 Conjugated Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Specificity	This antibody detects endogenous levels of UCHL1 and does not cross-react with related proteins.
Immunogen Description	Purified recombinant human UCHL1. protein fragments expressed in E.coli.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Neuron cytoplasmic protein 9.5; PGP 9.5; Ubiquitin thioesterase L1
Accession No.	Swiss-Prot#: P09936NCBI Gene ID:7345
Uniprot	P09936
GeneID	7345;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

## Background

---

Ubiquitin-protein hydrolase involved both in the processing of ubiquitin precursors and of ubiquitinated proteins. This enzyme is a thiol protease that recognizes and hydrolyzes a peptide bond at the C-terminal glycine of ubiquitin. Also binds to free monoubiquitin and may prevent its degradation in lysosomes. The homodimer may have ATP-independent ubiquitin ligase activity

---

Note: This product is for in vitro research use only