

Noggin Conjugated Monoclonal Antibody

Catalog No: #C27203



Package Size: #C27203-AF350 100ul #C27203-AF405 100ul #C27203-AF488 100ul
 #C27203-AF555 100ul #C27203-AF594 100ul #C27203-AF647 100ul
 #C27203-AF680 100ul #C27203-AF750 100ul #C27203-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	Noggin Conjugated Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Specificity	Transfected Only
Immunogen Description	Purified recombinant human NOG protein fragments expressed in E.coli.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	NOG; NOGG_HUMAN; Noggin; SYM 1; SYM1; Symphalangism 1 (proximal); Synostoses (multiple) syndrome 1; SYNS 1; SYNS1.
Accession No.	Swiss-Prot#: Q13253NCBI Gene ID:9241
Uniprot	Q13253
GeneID	9241;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

Essential for cartilage morphogenesis and joint formation. Inhibitor of bone morphogenetic proteins (BMP) signaling which is required for growth and patterning of the neural tube and somite.

Note: This product is for in vitro research use only