

Mouse anti-Human CD45, Purified mAbConjugated Antibody

Catalog No: #C28141

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C28141-AF350 100ul #C28141-AF405 100ul #C28141-AF488 100ul

#C28141-AF555 100ul #C28141-AF594 100ul #C28141-AF647 100ul

#C28141-AF680 100ul #C28141-AF750 100ul #C28141-Biotin 100ul

Description

Product Name	Mouse anti-Human CD45, Purified mAbConjugated Antibody
Host Species	Mouse
Clonality	Monoclonal
Species Reactivity	Hu
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

4AA9 reacts with CD45, a 180-220 kDa leukocyte common antigen (LCA). CD45 antigen is expressed at high levels on all hematopoietic cells including T and B lymphocytes, monocytes, granulocytes, NK cells and dendritic cells, but is not expressed on non-hematopoietic cells. CD45 antibody has also been reported to react weakly with mature blood erythrocytes and platelets. CD45 is a protein tyrosine phosphatase receptor that is critically important for T and B cell antigen receptor-mediated activation.

*Target of immunosuppressive antibody treatment.

Note: This product is for in vitro research use only