

CREB3L3 Antibody

Catalog No: #44850

Package Size: #44850-1 50ul #44850-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CREB3L3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	CREB3L3antibody detects endogenous levels of total?CREB3L3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CREB3L3
Target Name	CREB3L3
Other Names	cAMP responsive element binding protein 3 like 3 antibody cAMP responsive element-binding protein 3 like protein 3 antibody CREB H antibody CREB/ATF family transcription factor antibody CREBH antibody Cyclic AMP responsive element binding protein 3 like protein 3 antibody HYST1481 antibody MGC126553
Accession No.	Swiss-Prot#: Q68CJ9NCBI Gene ID:84699
Uniprot	Q68CJ9
GeneID	84699;
Calculated MW	49 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Note: This product is for in vitro research use only