

MIB2 Antibody

Catalog No: #44851

Package Size: #44851-1 50ul #44851-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MIB2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	MIB2antibody detects endogenous levels of total?MIB2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human MIB2
Target Name	MIB2
Other Names	MIB2
Accession No.	Swiss-Prot#: Q96AX9NCBI Gene ID:142678
Uniprot	Q96AX9
GeneID	142678;
Calculated MW	110 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Images



Western Blot analysis of HeLa lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only